



# Medical IT Trolley AMT200



**Dutch Design**  
**Medical IT**  
**of Tomorrow**

[www.athena-medical.com](http://www.athena-medical.com)

+31 85 876 9225  
[info@athena-medical.com](mailto:info@athena-medical.com)

# Medical IT Trolley

## AMT200



Extra strong internal Battery.  
Different drawer configurations available.

|                        |   |
|------------------------|---|
| Power system           | 220V AC input, 12V DC output                                      |
| Power status indicator | LED green-orange-red  |
| Monitor arm            | Gas spring  |
| Monitor weight         | 5 - 10kg  |
| VESA mount             | 75 x 75mm & 100 x 100mm   |
| Lift mechanism         | Motorized (DC 12V/Max. 7.2Amp/Push 500N)                          |
| Height adjustable      | Desktop 80cm to 110cm<br>Monitor 112cm to 153cm                   |
| Main material          | Aluminium alloy/Stainless steel/MDF with medic grade coating      |
| Mounting option        | DIN rails left and right (Option: DIN rails backsite)             |
| Work surface           | 510mm (L) x 580mm(W)  |
| Keyboard tray          | 260mm (D) x 536mm(W)  |
| Transport handle       | Front   |
| Wheels                 | 4 castors, 2 lockable castors in front<br>castor dimension: 125mm |
| Battery                | LiFePO 40Ah   |

|                      |  |
|----------------------|--|
| Battery compartment  | 180mm(L) x 290mm(W) x 159mm(H)   |
| Pre-assembled cables | 3 x 2.0 USB cables, 1 x LAN cable (RJ45 8P8C M to F)                   |
| Power cable          | 3 meter retractable power cable (220-240V/50-60Hz/6-8A)                |
| Compliance           | CE, RoHS, NEN/EN/IEC 3140, NEN/EN/IEC 62353, Class 2 double insulation |
| Security             | Integrated RFID reader (option)  |
| Warranty             | Mechanical parts: 2 years<br>Power system/electrical parts: 2 year     |



Acropolis Systems BV  
Ohmstraat 8  
3861 NB Nijkerk  
The Netherlands

+31 85 876 9225  
info@athena-medical.com  
www.athena-medical.com



## Medical IT Trolley **AMT200**

### **Ergonomic and Configurable Mobile Workstation**

The Athena AMT200 is an ergonomic and highly configurable Medical IT Trolley designed for demanding healthcare and high-standard environments. Available with both motorized and manually operated height adjustment, it supports comfortable seated and standing workflows while adapting to different ergonomic preferences.

Built for daily intensive use, the AMT200 combines stable mobility, robust construction, and long availability through an internal high-capacity battery. Its fully modular design allows configuration with drawers, additional work surfaces, device holders, and a medical-grade 24-inch All-in-One PC—making it a flexible and future-proof solution for clinical and medication-related workflows.



Acropolis Systems BV  
Ohmstraat 8  
3861 NB Nijkerk  
The Netherlands

+31 85 876 9225  
info@athena-medical.com  
www.athena-medical.com