



# Medical Trolley AMT500



**Dutch Design**  
**Medical IT**  
**of Tomorrow**

[www.athena-medical.com](http://www.athena-medical.com)

+31 85 876 9225  
[info@athena-medical.com](mailto:info@athena-medical.com)

# Medical Trolley AMT500



Power system	220V AC input, 12V DC output
Power status indicator	LED
Monitor arm	Gas spring
Monitor weight	5 - 10kg
VESA mount	75 x 75mm & 100 x 100mm
Lift mechanism	Gas spring (with foot pedal)
Height adjustable	Desktop 80cm to 110cm Monitor 112cm to 153cm
Main material	Aluminium alloy/Stainless steel/ABS plastic
Mounting option	DIN pole mounting for drawers and device holders
Work surface	450mm(L) x 500mm(W)
Keyboard tray	210mm(D) x 450mm(W) plus mouse table left/right
Transport handle	Front
Wheels	4 Castors, 2 lockable castors in front. Castor dimension: 125mm
Battery	Lithium Ion 15,6Ah / 230.9Wh

Battery compartment	Under table
Pre-assembled cables	3 x 2.0 USB cables 1 x LAN cable (RJ45 8P8C M to F)
Power cable	2 meter power cable (220-240V/50-60Hz/6-8A)
Compliance	CE, RoHS, NEN/EN/IEC 3140, NEN/EN/IEC 60335
Security (option)	Integrated RFID reader
Warranty	Mechanical parts: 2 years Power system/electrical parts: 2 year



Acropolis Systems BV  
Ohmstraat 8  
3861 NB Nijkerk  
The Netherlands

+31 85 876 9225  
info@athena-medical.com  
www.athena-medical.com



## Medical Trolley **AMT500**

### **Cost-effective and Compact Mobile Workstation**

The Athena AMT500 is an ergonomic and cost-effective Medical IT Trolley designed for environments where simplicity, mobility, and reliability are key. Featuring manual, hand-operated height adjustment, the AMT500 enables comfortable daily use with a compact, easy-to-handle design that supports ergonomic working positions.

Its modular structure allows integration of monitors, drawers, device holders, and accessories, while the internal battery supports mobile workflows throughout the day. The AMT500 is an ideal choice for departments seeking a financially attractive, ergonomic IT trolley that delivers dependable performance without the need for motorized height adjustment.



Acropolis Systems BV  
Ohmstraat 8  
3861 NB Nijkerk  
The Netherlands

+31 85 876 9225  
info@athena-medical.com  
www.athena-medical.com